
Nanoimager

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Extracellular vesicles (EVs)
细胞外囊泡

Virus research
Vaccine development
病毒研究、疫苗开发

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

20nm

RESOLUTION

Our combination of imaging modalities and advanced software allows single molecules to be observed and tracked with high sensitivity.

AUTOFOCUS

IN REALTIME

Rapid and precise stabilization of the microscope focus, with locking technology to minimize Z-drift.

CORRECTION

Post acquisition drift correction with nanometer precision

**TEMPERATURE
CONTROL**

**UP TO
42°**

Precise control and measurement of temperature, with heating elements to keep the entire microscope unit at 37°C.

Ideal for live cell imaging, facilitating single-particle tracking.

**POWERFUL
LASERS**

Single molecule localization microscopy in multiple colors and 2 channels.

Up to 1W laser power at the source.
8-12 kW/cm² power density at the sample plane.

**MICRO
FLUIDICS
COMPATIBLE**

Powerful solutions to control sample conditions.

DNA paint and RNA hybridization.

Up to 16 channels.

Automation of sample preparation.

Application

CASE STUDY 1

Enabling super-resolution microscopy imaging in biosafety cabinets

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
CASE STUDY 2


Single-molecule imaging for neuroscience

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Rat primary cortical neurons were cultured for 14 days, allowing synapse formation on a dish. Neuronal synapses were visualized in 3 color dSTORM imaging using two presynaptic proteins Bassoon (in magenta) and VGluT1 (in green); and the postsynaptic protein Homer1 (in orange).


Visualizing EV dynamics and interactions in living cells

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
CASE STUDY 5

Understanding the effects of viral infection

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
Rhinovirus (viral protein 2C-AF647, purple) infects HeLa cells and disrupts the Golgi (GM130-AF488, yellow). Acquired with Drs A Mousnier and G Schroeder, Queen's University Belfast.

CASE STUDY 6

How does the virus exploit the host cell machinery?

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dSTORM image of actin stained with PhalloidinAF647 and antibodies to nucleoprotein-AF555 in Ebola virus-like particle transfected cells.


- A) infected actin
- B) healthy cell
- C) virus protein nucleoprotein (red) colocalizes with disrupted actin network (blue)

CASE STUDY 7

Understanding pathogens with super-resolution microscopy

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Single-particle tracking of T3SS-HALO stained with JF549 in live bacteria.

CASE STUDY 8

Capturing organelle contact sites assembly to unveil disease onset in live cells

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▽ ER and MAMs (Mitochondria-associated ER membranes) protein PDZD8 (HaloTag®-PA, red) were imaged 'walking' along the ER in fibroblast cells in close vicinity of mitochondria (mitochondrial protein fused to GFP, green).

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△ Tracking of the MAM protein PDZD8 in fibroblast cells. (A) Individual protein traces were colored according to their diffusion coefficients. Panel on the right shows a selected trace overlaid with a widefield fluorescence signal from mitochondria (green).

Front cover: Atlastin1 protein walking along the ER in fibroblast cells. Sample provided by Dr. Christopher Obara, postdoc in Jennifer Lippincott-Schwartz lab, Jonelia Farm Applications page: CELL PHENOTYPING: dSTORM image of microtubules (green) and lamins (pink) *; NEUROSCIENCE: neuron labeled with AF568; BACTERIAL RESEARCH: single particle tracking of T3SS-HALO stained with JF549 in live bacteria, Dr. A. Diepold, Max Plank Marburg; IMMUNOLOGY: dSTORM of BCR-AF647 in B-cell, M-type, Dr. Bzostowski, NIH, Bethesda; VIRUS RESEARCH: dSTORM image of Marburg virus nucleoprotein labeled with AF647, prof. S. Becker, Marburg University; EXTRACELLULAR VESICLES: dSTORM image of an EV isolated from human keratinocyte culture media, CD63 (blue), CD81 (yellow), membrane surface (magenta) *; * Samples by ONI Applications Development Team

Back cover: dSTORM image of mitochondria in Cos7 cells stained with AF647 (Tomato20) and DyLight550 (Hsp60). Sample by ONI Imaging Team